

Human ECH1 recombinant protein



Catalog Number: 518320

General Information

Gene Name Synonym

Delta(3,5)-Delta(2,4)-dienoyl-CoA isomerase, mitochondrial, ECH1,, REF: C1008

Protein Construction

A DNA sequence encoding the Human ECH1 Thr34-Leu328 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris, 100mM NaCl, 10% Glycerol, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human ECH1 recombinant protein