

Human CXCL8/IL8 recombinant protein



Catalog Number: 518355

General Information

Gene Name Synonym

Interleukin-8, IL-8, C-X-C Motif Chemokine 8, Emotakin, Granulocyte Chemotactic Protein 1, GCP-1, Monocyte-Derived Neutrophil Chemotactic Factor, MDNCF, Monocyte-Derived Neutrophil-Activating Peptide, MONAP, Neutrophil-Activating Protein 1, NAP-1, Protein 3-10C, T-Cell Chemotactic Factor, IL8, CXCL8

Protein Construction

A DNA sequence encoding the Human CXCL8/IL8 Ser28-Ser99 is expressed.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted

samples are stable at < -20°C for 3 months.

Formulation

Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Usage Guide


Storage

Lyophilized protein should be stored at <-20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE



Human CXCL8/IL8 recombinant protein