

# Human SUMF1/FGE/PSEC0152/UNQ3037/PRO9852 (His tag) recombinant protein



Catalog Number: 518417

## General Information

### Gene Name Synonym

Sulfatase-Modifying Factor 1, C-Alpha-Formylglycine-Generating Enzyme 1, SUMF1, FGE

### Protein Construction

A DNA sequence encoding the Human SUMF1/FGE/PSEC0152/UNQ3037/PRO9852 Ser34-Asp374 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

## Stability

Store for up to 6 months at -20°C

## Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 2mM CaCl<sub>2</sub>, 10% Glycerol, pH 7.5.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human

SUMF1/FGE/PSEC0152/UNQ3037/PRO9852 (His tag) recombinant protein