

Human CSNK1G2/CK1G2 (His tag) recombinant protein



Catalog Number: 518421

General Information

Gene Name Synonym

Casein kinase I isoform gamma-2, CK1G2, CSNK1G2

Protein Construction

A DNA sequence encoding the Human CSNK1G2/CK1G2 Met18-Lys415 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM Tris, 500mM NaCl, 10% Glycerol, 1mM DTT, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human CSNK1G2/CK1G2 (His tag) recombinant protein