

# Mouse Ppia recombinant protein



Catalog Number: 518437

## General Information

### Gene Name Synonym

Peptidyl-prolyl cis-trans isomerase A; PPIase A; Cyclophilin A; Cyclosporin A-binding protein; Rotamase A; SP18; PPIA; CYPA

### Protein Construction

A DNA sequence encoding the Mouse Ppia Met1-Leu164 is expressed.

### Organism

Mouse

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2  $\mu$ m filtered solution of PBS, 10% glycerol, pH7.4.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Mouse Ppia recombinant protein