

# Human DEFB104A/DEFB104/DEFB4;/DEFB104B recombinant protein



Catalog Number: 518457

## General Information

### Gene Name Synonym

Beta-Defensin 104, Beta-Defensin 4, BD-4, DEFB-4, hBD-4, Defensin Beta 104, DEFB104A, DEFB104, DEFB4, DEFB104B

### Protein Construction

A DNA sequence encoding the Human DEFB104A/DEFB104/DEFB4;/DEFB104B Glu23-Pro72 is expressed.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

### Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

## Formulation

Lyophilized from a 0.2 µm filtered solution of 20mM PB, 130mM NaCl, pH 7.4.

## Usage Guide


### Storage

Lyophilized protein should be stored at <-20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

### Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## SDS-PAGE



Human DEFB104A/DEFB104/DEFB4;/DEFB104B recombinant protein