

Human PI4KB/PIK4CB (His tag) recombinant protein



Catalog Number: 518482

General Information

Gene Name Synonym

Phosphatidylinositol 4-Kinase Beta, PI4K-Beta, PI4KBeta, PtdIns 4-Kinase Beta, NPIK, PI4K92, PI4KB, PIK4CB

Protein Construction

A DNA sequence encoding the Human PI4KB/PIK4CB Met1-Met801 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 25mM Tris.HCl,pH7.3,100mM glycine,10%glycerol.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human PI4KB/PIK4CB (His tag) recombinant protein