

Human DYNLL1/DLC1/DNCL1/DNCLC1/HDLC1 (His tag) recombinant protein



Catalog Number: 518484

General Information

Gene Name Synonym

Dynein Light Chain 1 Cytoplasmic, 8 kDa Dynein Light Chain, DLC8, Dynein Light Chain LC8-Type 1, Protein Inhibitor of Neuronal Nitric Oxide Synthase, PIN, DYNLL1, DLC1, DNCL1, DNCLC1, HDLC1

Protein Construction

A DNA sequence encoding the Human DYNLL1/DLC1/DNCL1/DNCLC1/HDLC1 Met1-Gly89 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted

samples are stable at $< -20^{\circ}\text{C}$ for 3 months.

Formulation

Lyophilized from a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 200mM NaCl, 1mM DTT, pH 8.0.

Usage Guide


Storage

Lyophilized protein should be stored at $< -20^{\circ}\text{C}$, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at $4-7^{\circ}\text{C}$ for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than $100\ \mu\text{g}/\text{ml}$. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE



Human DYNLL1/DLC1/DNCL1/DNCLC1/HDLC1 (His tag) recombinant protein