

Human PRCP/PCP (His tag) recombinant protein



Catalog Number: 518495

General Information

Gene Name Synonym

Lysosomal Pro-X Carboxypeptidase,
Angiotensinase C, Lysosomal Carboxypeptidase C,
Proline Carboxypeptidase,
Prolylcarboxypeptidase, PRCP, PRCP, PCP

Protein Construction

A DNA sequence encoding the Human PRCP/PCP
Leu22-His496 is expressed with a 6His tag at the
C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL
method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM
NaAc-HAc, 150mM NaCl, 10% Glycerol, pH4.5.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human PRCP/PCP (His tag) recombinant protein