

# Human PSPH (His tag) recombinant protein



Catalog Number: 518505

## General Information

### Gene Name Synonym

Phosphoserine Phosphatase, PSP, PSPase, L-3-Phosphoserine Phosphatase, O-Phosphoserine Phosphohydrolase, PSPH

### Protein Construction

A DNA sequence encoding the Human PSPH Met1-Glu225 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 50mM Tris, 250mM NaCl, 50mM Imidazole, pH 8.5 .

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human PSPH (His tag) recombinant protein