

# Human HPCAL1/BDR1 (His tag) recombinant protein



Catalog Number: 518507

## General Information

### Gene Name Synonym

Hippocalcin-Like Protein 1, Calcium-Binding Protein BDR-1, HLP2, Visinin-Like Protein 3, VILIP-3, HPCAL1, BDR1

### Protein Construction

A DNA sequence encoding the Human HPCAL1/BDR1 Met1-Phe193 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

## Stability

Store for up to 6 months at -20°C

## Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 200mM NaCl, 1mM DTT, 10% Glycerol, pH 8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human HPCAL1/BDR1 (His tag) recombinant protein