

Human QPRT (His tag) recombinant protein



Catalog Number: 518523

General Information

Gene Name Synonym

Nicotinate-Nucleotide Pyrophosphorylase
[Carboxylating], Quinolinate
Phosphoribosyltransferase [Decarboxylating],
QAPRTase, QPRTase, QPRT

Protein Construction

A DNA sequence encoding the Human QPRT
Met1-His297 is expressed with a 6His tag at the
N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL
method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM
TrisHCl, 150mM NaCl, pH 8.0 .

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human QPRT (His tag) recombinant protein