

Human TRIM5/RNF88 (His tag) recombinant protein



Catalog Number: 518527

General Information

Gene Name Synonym

Tripartite motif-containing protein 5/RING finger protein 88/TRIM5/RNF88

Protein Construction

A DNA sequence encoding the Human TRIM5/RNF88 Met1-Gln248 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human TRIM5/RNF88 (His tag) recombinant protein