

# Mouse St6galnac2/Siat7/Siat7b (His tag) recombinant protein



Catalog Number: 518536

## General Information

### Gene Name Synonym

Gal-beta-1,3-GalNAc alpha-2,6-sialyltransferase, GalNAc alpha-2,6-sialyltransferase II, ST6GalNAc II, Sialyltransferase 7B, SIAT7-B

### Protein Construction

A DNA sequence encoding the Mouse St6galnac2/Siat7/Siat7b Ser28-Arg373 is expressed with a 6His tag at the C-terminus.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL

method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 50mM Tris, 150mM NaCl, pH7.5.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Mouse St6galnac2/Siat7/Siat7b (His tag)  
recombinant protein