

# Human CA7 (His tag) recombinant protein



Catalog Number: 518538

## General Information

### Gene Name Synonym

Carbonic Anhydrase 7, Carbonate Dehydratase VII, Carbonic Anhydrase VII, CA-VII, CA7

### Protein Construction

A DNA sequence encoding the Human CA7 Met1-Ala264 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL

method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0 .

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human CA7 (His tag) recombinant protein