

# Human SCLY/SCL (His tag) recombinant protein



Catalog Number: 518540

## General Information

### Gene Name Synonym

Selenocysteine Lyase, hSCL, SCLY, SCL

### Protein Construction

A DNA sequence encoding the Human SCLY/SCL Met1-Ala445 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

## Stability

Store for up to 6 months at -20°C

## Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human SCLY/SCL (His tag) recombinant protein