

# Human SULT2A1/HST/STD (His tag) recombinant protein



Catalog Number: 518546

## General Information

### Gene Name Synonym

Bile Salt Sulfotransferase, Dehydroepiandrosterone Sulfotransferase, DHEA-ST, Hydroxysteroid Sulfotransferase, HST, ST2, ST2A3, Sulfotransferase 2A1, ST2A1, SULT2A1, HST, STD

### Protein Construction

A DNA sequence encoding the Human SULT2A1/HST/STD Ser2-Glu285 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM Tris, 100mM NaCl, pH 8.0.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human SULT2A1/HST/STD (His tag) recombinant protein