

# Human PPIA/CYPA recombinant protein



Catalog Number: 518548

## General Information

### Gene Name Synonym

Peptidyl-prolyl cis-trans isomerase A; PPIase A; Cyclophilin A; Cyclosporin A-binding protein; Rotamase A; SP18; PPIA; CYPA

### Protein Construction

A DNA sequence encoding the Human PPIA/CYPA Met1-Glu165 is expressed.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL

method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of PBS, 10% glycerol, pH7.4.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human PPIA/CYPA recombinant protein