

Human MIF/GLIF/MMIF recombinant protein



Catalog Number: 518559

General Information

Gene Name Synonym

Macrophage migration inhibitory factor, MIF, MMIF, Glycosylation-inhibiting factor, GLIF, L-dopachrome tautomerase, Phenylpyruvate tautomerase

Protein Construction

A DNA sequence encoding the Human MIF/GLIF/MMIF Met1-Ala115 is expressed.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human MIF/GLIF/MMIF recombinant protein