

## General Information

### Gene Name Synonym

Macrophage migration inhibitory factor, MIF, MMIF, Glycosylation-inhibiting factor, GLIF, L-dopachrome tautomerase, Phenylpyruvate tautomerase

### Protein Construction

A DNA sequence encoding the Human MIF/GLIF/MMIF Met1-Ala115 is expressed.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human MIF/GLIF/MMIF recombinant protein