

Human RCN3/UNQ239/PRO272 (His tag) recombinant protein



Catalog Number: 518581

General Information

Gene Name Synonym

Reticulocalbin-3, EF-Hand Calcium-Binding Protein RLP49, RCN3

Protein Construction

A DNA sequence encoding the Human RCN3/UNQ239/PRO272 Lys21-Leu328 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 150mM NaCl, 1mM DTT, 10% Glycerol, pH7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human RCN3/UNQ239/PRO272 (His tag)
recombinant protein