

# Human PRDX3/AOP1 recombinant protein



Catalog Number: 518591

## General Information

### Gene Name Synonym

Thioredoxin-Dependent Peroxide Reductase Mitochondrial, Antioxidant Protein 1, AOP-1, HBC189, Peroxiredoxin III, Prx-III, Peroxiredoxin-3, Protein MER5 homolog, PRDX3, AOP1

### Protein Construction

A DNA sequence encoding the Human PRDX3/AOP1 Pro63-Gln256 is expressed.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM TrisHCl, pH8.0.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human PRDX3/AOP1 recombinant protein