

Human DCK (His tag) recombinant protein



Catalog Number: 518658

General Information

Gene Name Synonym

Deoxycytidine Kinase, dCK, DCK

Protein Construction

A DNA sequence encoding the Human DCK Met1-Leu260 is expressed with a 6His, T7 tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl pH7.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human DCK (His tag) recombinant protein