

# Human ESD (His tag) recombinant protein



Catalog Number: 518665

## General Information

### Gene Name Synonym

S-Formylglutathione Hydrolase, FGH, Esterase D, Methylumbelliferyl-Acetate Deacetylase, ESD

### Protein Construction

A DNA sequence encoding the Human ESD Met1-Ala282 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2  $\mu$ m filtered solution of 20mM TrisHCl, 10% Glycerol, pH 8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human ESD (His tag) recombinant protein