

# Human CD86/CD28LG2 (His tag) recombinant protein



Catalog Number: 518694

## General Information

### Gene Name Synonym:

T-Lymphocyte Activation Antigen CD86, Activation B7-2 Antigen, B70, BU63, CTLA-4 Counter-Receptor B7.2, FUN-1, CD86, CD28LG2

### Protein Construction:

A DNA sequence encoding the Human CD86/CD28LG2 Ala24-Pro247 is expressed with a 6His tag at the C-terminus.

**Source:** Human

**Expression Host:** Human Cells

## QC Testing

**Activity:**

**Purity:**

>95% as determined by SDS-PAGE.

**Endotoxin:**

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

**Stability:**

Store for up to 12 months at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  as lyophilized powder. Aliquots of reconstituted samples are stable at  $< -20^{\circ}\text{C}$  for 3 months.

**Predicted N terminal:**

**Molecular Mass:**

## Formulation:

Lyophilized from a  $0.2\ \mu\text{m}$  filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

## Usage Guide


### Storage:

Lyophilized protein should be stored at  $< -20^{\circ}\text{C}$ , though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at  $4-7^{\circ}\text{C}$  for 2-7 days.

### Reconstitution:

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than  $100\ \mu\text{g}/\text{ml}$ . Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

### SDS-PAGE:



Human CD86/CD28LG2 (His tag) recombinant protein