

Human LYG2/LYGH (His tag) recombinant protein



Catalog Number: 518736

General Information

Gene Name Synonym

Lysozyme G-Like Protein 2, LYG2, LYGH

Protein Construction

A DNA sequence encoding the Human LYG2/LYGH Ser20-Phe212 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, 10% Glycerol, pH 7.2.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human LYG2/LYGH (His tag) recombinant protein