

# Mouse Nt5e/Nt5/Nte (His tag) recombinant protein



Catalog Number: 518802

## General Information

### Gene Name Synonym

5'-nucleotidase, Ecto-5'-nucleotidase, CD73, 5'-NT

### Protein Construction

A DNA sequence encoding the Mouse Nt5e/Nt5/Nte Trp29-Phe550 is expressed with a 6His tag at the C-terminus.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

## Stability

Store for up to 6 months at -20°C

## Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 120mM NaCl, 4mM CaCl<sub>2</sub>, 20% Glycerol, pH 7.5.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

## SDS-PAGE



Mouse Nt5e/Nt5/Nte (His tag) recombinant protein