

# Human S100A8/S100A9 Heterodimer (His tag) recombinant protein



Catalog Number: 518826

## General Information

### Gene Name Synonym

S100A8/S100A9 Heterodimer

### Protein Construction

A DNA sequence encoding the Human S100A8/S100A9 Heterodimer Met1-Glu93(S100A8)&Thr2-Pro114(S100A9) is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

## Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

## Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM HEPES, 10% Glycerol, 150mM NaCl, 2.5mM EDTA, pH 7.4.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

## SDS-PAGE



Human S100A8/S100A9 Heterodimer (His tag)  
recombinant protein