

Human S100A8/S100A9 Heterodimer (His tag) recombinant protein



Catalog Number: 518826

General Information

Gene Name Synonym

S100A8/S100A9 Heterodimer

Protein Construction

A DNA sequence encoding the Human S100A8/S100A9 Heterodimer Met1-Glu93(S100A8)&Thr2-Pro114(S100A9) is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM HEPES, 10% Glycerol, 150mM NaCl, 2.5mM EDTA, pH 7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human S100A8/S100A9 Heterodimer (His tag)
recombinant protein