

Human HAMP/HEPC/LEAP1/UNQ487/PRO1003 (GST tag) recombinant protein



Catalog Number: 518827

General Information

Gene Name Synonym

Liver-expressed antimicrobial peptide 1, Putative liver tumor regressor, HEPC, LEAP1

Protein Construction

A DNA sequence encoding the Human HAMP/HEPC/LEAP1/UNQ487/PRO1003 Asp60-Thr84 is expressed with a GST tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of PBS□50%glycerol□pH7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human HAMP/HEPC/LEAP1/UNQ487/PRO1003 (GST tag) recombinant protein