

Human BSG/UNQ6505/PRO21383 (His tag) recombinant protein



Catalog Number: 518841

General Information

Gene Name Synonym

Basigin, 5F7, Collagenase Stimulatory Factor, Extracellular Matrix Metalloproteinase Inducer, EMMPRIN, Leukocyte Activation Antigen M6, OK Blood Group Antigen, Tumor Cell-Derived Collagenase Stimulatory Factor, TCSF, CD147, BSG

Protein Construction

A DNA sequence encoding the Human BSG/UNQ6505/PRO21383 Ala22-His205 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted

samples are stable at $< -20^{\circ}\text{C}$ for 3 months.

Formulation

Lyophilized from a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Usage Guide


Storage

Lyophilized protein should be stored at $< -20^{\circ}\text{C}$, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at $4-7^{\circ}\text{C}$ for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than $100\ \mu\text{g}/\text{ml}$. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE



Human BSG/UNQ6505/PRO21383 (His tag)
recombinant protein