

Human GOT1 (His tag) recombinant protein



Catalog Number: 518848

General Information

Gene Name Synonym

Aspartate Aminotransferase Cytoplasmic, Glutamate Oxaloacetate Transaminase 1, Transaminase A, GOT1

Protein Construction

A DNA sequence encoding the Human GOT1 Ala2-Leu414 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 100mM NaCl, 2mM DTT, 20% Glycerol, pH 7.5 .

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human GOT1 (His tag) recombinant protein