

# Mouse Mmp9/Clg4b recombinant protein



Catalog Number: 518889

## General Information

### Gene Name Synonym

Matrix metalloproteinase-9; MMP-9; 92 kDa gelatinase; 92 kDa type IV collagenase; Gelatinase B; GELB

### Protein Construction

A DNA sequence encoding the Mouse Mmp9/Clg4b Ala20-Pro730 is expressed with a 10His tag at the C-terminus.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH7.5.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Mouse Mmp9/Clg4b recombinant protein