

# Human KLK11/PRSS20/TLSP/UNQ649/PRO1279 (His tag) recombinant protein



Catalog Number: 518911

## General Information

### Gene Name Synonym

Kallikrein-11, hK11, Hippostasin, Serine Protease 20, Trypsin-Like Protease, KLK11, PRSS20, TLSP

### Protein Construction

A DNA sequence encoding the Human KLK11/PRSS20/TLSP/UNQ649/PRO1279 Glu19-Asn250 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 2mM CaCl<sub>2</sub>, pH 8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human KLK11/PRSS20/TLSP/UNQ649/PRO1279 (His tag) recombinant protein