

Human CD40LG/CD40L/TNFSF5/TRAP recombinant protein



Catalog Number: 518954

General Information

Gene Name Synonym

CD40 Ligand, CD40-L, T-Cell Antigen Gp39, TNF-Related Activation Protein, TRAP, Tumor Necrosis Factor Ligand Superfamily Member 5, CD154, CD40LG, CD40L, TNFSF5, TRAP

Protein Construction

A DNA sequence encoding the Human CD40LG/CD40L/TNFSF5/TRAP Met113-Leu261 is expressed.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Formulation

Lyophilized from a 0.2 µm filtered solution of 20mM PB, 200mM NaCl, 0.1mM EDTA, pH 7.0.

Usage Guide


Storage

Lyophilized protein should be stored at <-20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE



Human CD40LG/CD40L/TNFSF5/TRAP recombinant protein