

Human KLK2 (His tag) recombinant protein



Catalog Number: 519037

General Information

Gene Name Synonym

Kallikrein-2; Glandular Kallikrein-1; hGK-1; Tissue Kallikrein-2; KLK2

Protein Construction

A DNA sequence encoding the Human KLK2 Pro19-Pro261 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM Citrate, 150mM NaCl, pH 3.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human KLK2 (His tag) recombinant protein