

# Human GALNT7 (His tag) recombinant protein



Catalog Number: 519041

## General Information

### Gene Name Synonym

PolypeptideGalNAc-transferase7,Protein-UDP acetylgalactosaminyltransferase7,UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase 7, GALNT7.

### Protein Construction

A DNA sequence encoding the Human GALNT7 Pro30-Val657 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2  $\mu$ m filtered solution of 50mM TrisHCl,10mM reduced Glutathione,PH8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human GALNT7 (His tag) recombinant protein