

Human GALNT7 (His tag) recombinant protein



Catalog Number: 519041

General Information

Gene Name Synonym

PolypeptideGalNAc-transferase7,Protein-UDP acetylgalactosaminyltransferase7,UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase 7, GALNT7.

Protein Construction

A DNA sequence encoding the Human GALNT7 Pro30-Val657 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 50mM TrisHCl, 10mM reduced Glutathione, PH8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human GALNT7 (His tag) recombinant protein