

Human LTA4H/LTA4 (His tag) recombinant protein



Catalog Number: 519051

General Information

Gene Name Synonym

Leukotriene A-4 hydrolase, LTA-4 hydrolase, EC=3.3.2.6, Leukotriene A(4) hydrolase, LTA4

Protein Construction

A DNA sequence encoding the Human LTA4H/LTA4 Met1-Asp611 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM MES, 2mM EDTA, 20% Glycerol, pH 5.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human LTA4H/LTA4 (His tag) recombinant protein