

Human UBE2K/HIP2/LIG (GST tag) recombinant protein



Catalog Number: 519102

General Information

Gene Name Synonym

Ubiquitin-Conjugating Enzyme E2 K, Huntingtin-Interacting Protein 2, HIP-2, Ubiquitin Carrier Protein, Ubiquitin-Conjugating Enzyme E2-25 kDa, Ubiquitin-Conjugating Enzyme E2(25K), Ubiquitin-Conjugating Enzyme E2-25K, Ubiquitin-Protein Ligase, UBE2K, HIP2, LIG

Protein Construction

A DNA sequence encoding the Human UBE2K/HIP2/LIG Met1-Asn200 is expressed with a GST tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human UBE2K/HIP2/LIG (GST tag) recombinant protein