

Human AGR2/AG2/UNQ515/PRO1030 (His tag) recombinant protein



Catalog Number: 519180

General Information

Gene Name Synonym

Anterior Gradient Protein 2 Homolog, AG-2, hAG-2, HPC8, Secreted Cement Gland Protein XAG-2 Homolog, AGR2, AG2

Protein Construction

A DNA sequence encoding the Human AGR2/AG2/UNQ515/PRO1030 Arg21-Leu175 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris, 200mM NaCl, 10% Glycerol, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human AGR2/AG2/UNQ515/PRO1030 (His tag)
recombinant protein