

Human CCL23/MIP3/MPIF1/SCYA23 recombinant protein



Catalog Number: 519189

General Information

Gene Name Synonym

C-C Motif Chemokine 23, CK-Beta-8, CKB-8, Macrophage Inflammatory Protein 3, MIP-3, Myeloid Progenitor Inhibitory Factor 1, MPIF-1, Small-Inducible Cytokine A23, CCL23, MIP3, MPIF1, SCYA23

Protein Construction

A DNA sequence encoding the Human CCL23/MIP3/MPIF1/SCYA23 Arg22-Asn120 is expressed.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted samples are stable at $< -20^{\circ}\text{C}$ for 3 months.

Formulation

Lyophilized from a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 250mM NaCl, pH 7.2.

Usage Guide


Storage

Lyophilized protein should be stored at $< -20^{\circ}\text{C}$, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at $4-7^{\circ}\text{C}$ for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than $100\ \mu\text{g}/\text{ml}$. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE



Human CCL23/MIP3/MPIF1/SCYA23 recombinant protein