

Human UBE2I/UBC9/UBCE9 (GST tag) recombinant protein



Catalog Number: 519231

General Information

Gene Name Synonym

SUMO-Conjugating Enzyme UBC9, SUMO-Protein Ligase, Ubiquitin Carrier Protein 9, Ubiquitin Carrier Protein I, Ubiquitin-Conjugating Enzyme E2 I, Ubiquitin-Protein Ligase I, p18, UBE2I, UBC9, UBCE9

Protein Construction

A DNA sequence encoding the Human UBE2I/UBC9/UBCE9 Met1-Ser158 is expressed with a GST tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 50mM HEPES, 150mM NaCl, pH 7.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human UBE2I/UBC9/UBCE9 (GST tag)
recombinant protein