

Human LRG1/LRG (Fc & His tag) recombinant protein



Catalog Number: 519247

General Information

Gene Name Synonym

Leucine-rich alpha-2-glycoprotein, LRG and LRG1

Protein Construction

A DNA sequence encoding the Human LRG1/LRG Val36-Gln347 is expressed with a Fc, 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris, 200mM NaCl, 10% Glycerol pH8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human LRG1/LRG (Fc & His tag) recombinant protein