

# Human SDCBP/MDA9/SYCL (His tag) recombinant protein



Catalog Number: 519250

## General Information

### Gene Name Synonym

Syntenin-1, Melanoma differentiation-associated protein 9, Pro-TGF-alpha cytoplasmic domain-interacting protein 18, Scaffold protein Pbp1, Syndecan-binding protein 1, SDCBP, MDA9, SYCL,

### Protein Construction

A DNA sequence encoding the Human SDCBP/MDA9/SYCL Ser2-Val298 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM PB, 150mM NaCl, pH7.4.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human SDCBP/MDA9/SYCL (His tag) recombinant protein