

# Mouse Pla2g1b/Pla2 (His tag) recombinant protein



Catalog Number: 519259

## General Information

### Gene Name Synonym

Phospholipase A2; Group IB phospholipase A2; PLA2-Ib; Phosphatidylcholine 2-acylhydrolase 1B; Pla2g1b; Pla2

### Protein Construction

A DNA sequence encoding the Mouse Pla2g1b/Pla2 Ala16-Cys146 is expressed with a 6His tag at the C-terminus.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM HEPES, 150mM NaCl, pH 7.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Mouse Pla2g1b/Pla2 (His tag) recombinant protein