

Human HTN3/HIS2 recombinant protein



Catalog Number: 519358

General Information

Gene Name Synonym

Histatin-3; Basic histidine-rich protein; Histidine-rich protein 3; PB; HTN3.

Protein Construction

A DNA sequence encoding the Human HTN3/HIS2 Asp20-Asn51 is expressed.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of PBS, pH 7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human HTN3/HIS2 recombinant protein