

Mouse Scarb2 (His tag) recombinant protein



Catalog Number: 519390

General Information

Gene Name Synonym

Lysosome membrane protein 2, 85 kDa lysosomal membrane sialoglycoprotein, LGP85, Lysosome membrane protein II, LIMP II, Scavenger receptor class B member 2, SCARB2

Protein Construction

A DNA sequence encoding the Mouse Scarb2 Arg27-Thr432 is expressed with a Fc tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Formulation

Lyophilized from a 0.2 µm filtered solution of 50mM Tris-Citrate, 0.3M NaCl, pH6.5.

Usage Guide


Storage

Lyophilized protein should be stored at <-20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE



Mouse Scarb2 (His tag) recombinant protein