

# Human UBE2B/RAD6B (His tag) recombinant protein



Catalog Number: 519401

## General Information

### Gene Name Synonym

Ubiquitin-Conjugating Enzyme E2 B, RAD6 Homolog B, HR6B, hHR6B, Ubiquitin Carrier Protein B, Ubiquitin-Conjugating Enzyme E2-17 kDa, Ubiquitin-Protein Ligase B, UBE2B, RAD6B

### Protein Construction

A DNA sequence encoding the Human UBE2B/RAD6B Met1-Ser152 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2  $\mu$ m filtered solution of 20mM PB, 150mM NaCl, pH7.2.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human UBE2B/RAD6B (His tag) recombinant protein