

# Human CADM3/IGSF4B/NECL1 (His tag) recombinant protein



Catalog Number: 519412

## General Information

### Gene Name Synonym

Cell Adhesion Molecule 3, Brain Immunoglobulin Receptor, Immunoglobulin Superfamily Member 4B, IgSF4B, Nectin-Like Protein 1, NECL-1, Synaptic Cell Adhesion Molecule 3, SynCAM3, TSLC1-Like Protein 1, TSLL1, CADM3, IGSF4B, NECL1, SYNCAM3, TSLL1

### Protein Construction

A DNA sequence encoding the Human CADM3/IGSF4B/NECL1 Asn25-His330 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

### Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted

samples are stable at < -20°C for 3 months.

### Formulation

Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

## Usage Guide


### Storage

Lyophilized protein should be stored at <-20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

### Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## SDS-PAGE



Human CADM3/IGSF4B/NECL1 (His tag)  
recombinant protein